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
*The Squidstat Solo and Boosted Squidstat series don't use an external power supply.

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View All

Search

- Charge/Discharge 1 (CC-CV)
- Charge/Discharge 2 (CC-CV)_n-PotEIS
- Charge/Discharge 3 (CC-CV)_n-GalvEIS
- Charge/Discharge 4 (CP/CR)
- Chronoamperometry
- Chronopotentiometry
- Constant Power
- Constant Resistance
- Cyclic Voltammetry
- Differential Pulse Voltammetry
- Galvanostatic EIS
- Potentiostatic EIS
- GITT

Parameters

Quiet time: 0
(time spent at the starting potential)

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
Constant Voltage Stage

Rest Phase

CC-CV Charge/Discharge (Galvanostatic Cycling with Potential Limitation)

This experiment subjects a battery to constant current-constant voltage (CC-CV) charge/discharge cycles. For each charging and discharging phase, the battery undergoes a constant current (CC) stage and then a constant voltage (CV) stage. The CC stage ends when the Voltage limit is reached. The CV stage ends when the Current limit is reached. The device can be put at a rest phase after every charge or discharge phase. When the device is at rest, no current is flowing through.

Parameters

General Options

Starting phase: Charge first

Skip constant voltage steps? Yes No

Cycles: 1

Start a new data file for each cycle? Yes No

Create new file on charge and discharge of each cycle? Yes No

Charging Options

Current: 100 mA

Voltage limit: 4.2 V

End constant voltage step when [current] <: 10 mA

Constant Current (CC) sampling interval: 1 s

Constant Voltage (CV) sampling interval: 1 s

Constant Current (CC) time limit: 300 s

Constant Voltage (CV) time limit: 500 s

CC - CV time limit: Optional s

Discharging Options

Current: 100 mA

Voltage limit: 2.8 V

End constant voltage step when [current] <: 10 mA

Constant Current (CC) sampling interval: 1 s

Constant Voltage (CV) sampling interval: 1 s

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Run an Experiment	View Data	Multichannel Control	Build an Experiment	Manual Control	More Options
-------------------	-----------	----------------------	---------------------	----------------	--------------

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